Atellica IM¹

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3. Boeddighaus J, Twerenbold R, Nestelberger T, et al. Clinical validation of a novel high-sensitivity cardiac troponin I assay for early diagnosis of acute myocardial infarction. *Clin Chem.* 2018;64:1347-60.
4. Morrow DA, Cannon CP, Jesse RL, et al. National Academy of Clinical Biochemistry Laboratory Medicine Practice Guidelines: clinical characteristics and utilization of biochemical markers in acute coronary syndromes. *Clin Chem.* 2007;53:552-74.
5. Myocardial infarction redefined—a consensus document of The Joint European Society of Cardiology/American College of Cardiology Committee for the Redefinition of Myocardial Infarction. *Eur Heart J.* 2000;21:1502-13.
6. Apple FS, Collinson PO. Analytical characteristics of high-sensitivity cardiac troponin assays. *Clin Chem.* 2012;58:54-61.
7. Hamm CW, Bassand JP, Agewall S, et al. ESC Guidelines for the management of acute coronary syndromes in patients presenting without persistent ST-segment elevation: The Task Force for the Management of Acute Coronary Syndromes (ACS) in Patients Presenting without Persistent ST-segment Elevation of the European Society of Cardiology (ESC). *Eur Heart J.* 2011;32:2999-3054.
8. Thygesen K, Alpert JS, Jaffe AS, et al. Fourth universal definition of myocardial infarction. *Eur Heart J.* 2018;00:1-33.
9. Twerenbold R, Boeddighaus J, Mueller C. Update on high-sensitivity cardiac troponin in patients with suspected myocardial infarction. *Eur Heart J. Suppl* 2018;(Supplement G):G2-G10.
10. Keller T, Zeller T, Peetz D, et al. Sensitive troponin I assay in early diagnosis of acute myocardial infarction. *N Engl J Med.* 2009;361:868-77.
11. Reichlin T, Hochholzer W, Bassetti S, et al. Early diagnosis of myocardial infarction with sensitive cardiac troponin assays. *N Engl J Med.* 2009;361:858-67.
12. Wu AH, Lu QA, Todd J, Moeckls J, Wians F. Short- and long-term biological variation in cardiac troponin I measured with a high-sensitivity assay: implications for clinical practice. *Clin Chem.* 2009;55:52-8.
13. Reichlin T, Irfan A, Twerenbold R, et al. Utility of absolute and relative changes in cardiac troponin concentrations in the early diagnosis of acute myocardial infarction. *Circulation.* 2011;124:136-45.
14. Jaffe AS, Moeckel M, Giannitsis E, et al. In search for the Holy Grail: suggestions for studies to define delta changes to diagnose or exclude acute myocardial infarction: a position paper from the Study Group on Biomarkers of the Acute Cardiovascular Care Association. *European Heart Journal - Acute Cardiovascular Care.* 2014;3:313-6.
15. Jaffe AS, Babuin L, Apple FS. Biomarkers in acute cardiac disease: the present and the future. *J Am Coll Cardiol.* 2006;48:1-11.
16. Scharnhorst V, Krasznai K, van 't Veer M, Michels RH. Variation of cardiac troponin I and T measured with sensitive assays in emergency department patients with noncardiac chest pain. *Clin Chem.* 2012;58:1208-14.
17. Thygesen K, Mair J, Giannitsis E, et al. How to use high-sensitivity cardiac troponins in acute cardiac care. *Eur Heart J.* 2012;33:2252-7.
18. Wu AH, Bolger AF, Hollander JE. Growing pains with the use of high-sensitivity cardiac troponin assays. *J Am Coll Cardiol.* 2013;62:1250-1.
19. Clinical and Laboratory Standards Institute. CLSI document EP09-A3: Measurement procedure comparison and bias estimation using patient samples. Approved guideline-third edition. CLSI, Wayne, PA. 2013.
20. Roffi M, Patrono C, Collet JP, et al. 2015 ESC Guidelines for the management of acute coronary syndromes in patients presenting without persistent ST-segment elevation: Task Force for the Management of Acute Coronary Syndromes in Patients Presenting without Persistent ST-Segment Elevation of the European Society of Cardiology (ESC). *Eur Heart J.* 2016;37:267-315.
21. Stepinska J, et al. *Eur Heart J Acute Cardiovasc Care.* 2020 Feb;9(1):76-89.
22. Chapman AR, et al. *Heart* 2018;0:1–7. doi:10.1136/heartjnl-2018-314093.
23. Nowak RM, Christenson RH, Jacobsen G, et al. Performance of novel high-sensitivity cardiac troponin I assays for 0/1-hour and 0/2- to 3-hour evaluations for acute myocardial infarction: Results from the HIGH-US study. *Ann Emerg Med* 2020;76:1-13.

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Siemens Healthcare GmbH